

Lab 1.

BASIC SKILLS: DILUTIONS, MICROPIPETTES AND SPECTROPHOTOMETRY

There are three parts to this exercise:

- I. Concepts and Calculations for Dilutions and Solutions
- II. Use and Calibration of Micropipettes
- III. Spectrophotometric Analysis of Food Colors

I. CONCEPTS AND CALCULATIONS FOR SOLUTIONS

Preparing reagents and solutions is a never-ending task in most laboratories. This is a basic laboratory skill that often confuses people at first. Here we present the standard, general approach to computing dilutions and concentrations; the Dilution Factor Technique. It is a convenient way of computing dilutions at the bench. Work through this section BEFORE coming to lab.

A. Terminology and Concepts

Stock solution: concentrated solution which is being diluted

Working solution: diluted solution, ready to use

Diluent: the fluid used for diluting concentrate

Dilution Factor (DF):

The DF is a dimensionless number that unambiguously describes the "strength" of the dilution. It is equal to the volume of stock solution used (V_1), divided by the total volume of working solution produced (V_2). In turn, $V_2 = V_1 +$ the volume of diluent used.

The dilution factor also gives the relationship between solute concentration in the stock solution (C_1) and the working solution (C_2).

The precise relationships are given by:

$$\text{Dilution factor (DF)} = C_2 / C_1 = V_1 / V_2$$

This is sometimes given as:

$$C_1 \times V_1 = C_2 \times V_2$$

- Remember:
- $V_2 > V_1$
 - $C_2 < C_1$
 - $DF < 1.0$ (A "50-fold [50X]" dilution has $DF = 1/50 = 0.02$)
 - When you add 1.0 ml of stock to 4.0 ml of diluent, $DF = 1/5$, NOT $1/4$.

Serial dilutions

A sequential set of dilutions in which the stock for each dilution in the series is the working solution from previous dilution. In effect, except for the last dilution, each dilution is both a stock and a working solution.

The DF for the entire series as a whole is the product of the DF's of each individual dilution.

$$DF_{\text{total}} = DF_1 \times DF_2 \times DF_3 \text{ etc.}$$

Ordinarily, volumes are chosen to give DF values that are simple powers of ten. This makes it relatively simple to remember what you are doing and to perform the dilution calculations in your head.

The most commonly encountered standard dilutions are described in the table below:

Stock Volume (V1)	Diluent Volume	Working Volume (V2)	DF
1.0 ml	9.0 ml	10.0 ml	10^{-1}
0.1 ml	0.9 ml	1.0 ml	10^{-1}
0.1 ml	9.9 ml	10.0 ml	10^{-2}
0.01 ml	0.99 ml	1.0 ml	10^{-2}

B. EXAMPLES

1. Given a stock of 32% w/v NaCl, how would you make 500 ml of 5.6% w/v NaCl?

- $DF = C_2 / C_1 = 5.6/32 = 0.175$
- $V_1 = DF \times V_2 = 0.175 \times 500 \text{ ml} = 87.5 \text{ ml}$
- Add 87.5 ml NaCl stock to 412.5 ml water.

2. How would you make 100 ml of 32% w/v NaCl? (Note: A 1% w/v solution contains 1 g of solute in a volume of 100 ml).

- Weigh 32.0 g NaCl on a balance.
- Add enough water to give a final volume of 100 ml. (This will be less than 100 ml!)

3. You have a stock solution of 20X TAE Buffer. "20X" refers to the solute concentration of the stock relative to the final "working concentration". How would you prepare 1 liter of TAE Buffer at working concentration (1X) from the 20X TAE stock?

- $DF = C_2 / C_1 = 1X/20X = 0.05$
- $V_1 = DF \times V_2 = 0.05 \times 1,000 \text{ ml} = 50.0 \text{ ml}$
- Measure 50.0 ml of the 20X stock solution and add enough water to give a final volume of 1 liter. (The approximate volume of water will be 950 ml.)

4. How much of a 20 mg/ml stock solution of ampicillin would you need to add to 1 liter of bacterial growth medium to achieve a final concentration of 100 µg/ml ampicillin?

- $DF = C_2 / C_1 = 100 \text{ µg/ml} / 20,000 \text{ µg/ml} = 5 \times 10^{-3}$
- $V_1 = DF \times V_2 = 5 \times 10^{-3} \times 1,000 \text{ ml} = 5.0 \text{ ml}$
- Add 5.0 ml ampicillin stock to 1 liter medium.

5. You have 1.0 ml of bovine serum albumin (BSA) at a concentration of 10 mg/ml. You want 1.0 ml of BSA at a concentration of 10 ng/ml. Describe the procedure you would use to dilute your sample to 10 ng/ml.

Careful. This is not as straightforward as it appears. Remember that you cannot measure volumes less than 1 μl accurately; therefore a serial dilution is used.

- $DF = C_2 / C_1 = 10 \text{ ng/ml} / 10^7 \text{ ng/ml} = 10^{-6}$
- $V_1 = DF \times V_2 = 10^{-6} \times 1,000 \mu\text{l} = 10^{-3} \mu\text{l}$
- Add $10^{-3} \mu\text{l}$ stock to 1.0 ml water.
- The problem we see with this method is that there is no convenient way to measure $10^{-3} \mu\text{l}$ accurately.
- Instead, you should perform a series of standard dilutions. For example, could do a series of three consecutive 10^{-2} dilutions to reach a total dilution of 10^{-6} .

C. Exercises

1. Fill in the blank spaces with the appropriate value. See Appendix A.7 for the definitions of unit prefixes in the metric system.

10 mM = _____ M
1 μ M = _____ M
1 liter = _____ ml
0.1 ml = _____ μ l
10 mg = _____ g
10 μ g = _____ g

2. Describe how you would prepare 200 ml of a 0.25 M solution of NaOH (molecular weight = 40 g/mol).
3. "Physiological" saline is a solution of 0.85% (weight/volume) NaCl that is isotonic with tissue. Describe how you would prepare 1 liter of physiological saline. (Remember that 1% w/v is equivalent to 1 g/100 ml.)
4. Calculate how much glycerol you would need to make 200 ml of a 20% (v/v) glycerol solution. (Glycerol is a liquid.)
5. Describe how you would prepare 125 ml of a 0.8% w/v agarose in 1X TAE buffer starting with a 20X TAE stock solution and solid agarose powder.
6. If a 100X solution is 1 M, what is the concentration of a 1X solution?
7. The following problem is derived from a notorious article in a prestigious journal. [Davenas, E. *et al.* (1988) *Nature* 333, 816-818.].

Solutions containing a specific antibody cause a response called "degranulation" when the antibody solution is mixed with a suspension of human leukocytes (white blood cells). Degranulation is part of the normal immunologic response of the leukocytes to the presence of foreign antigens. In an assay for degranulation, 1.5 ml of the antibody solution is mixed with 30 μ l of leukocyte cell suspension and then the cells are examined using a microscopic to determine whether or not degranulation has occurred.


The article in *Nature* reported the following experimental observation based on the degranulation assay described above: A solution of antibody at a concentration of 2×10^{-9} M causes the degranulation response. When the 2×10^{-9} M antibody solution was diluted by a factor of $DF = 10^{-100}$, the degranulation response still occurred.

Why do you think this claim evoked expressions of astonishment worldwide? Evaluate the result using the ability to perform dilution series calculations and suggest why it was met with incredulity. (HINT: Avagadro's number is 6.02×10^{23} molecules/mole.)

II. USE AND CALIBRATION OF MICROPIPETTES

Take nothing for granted, least of all the equipment you rely on. Micropipettes are indispensable tools in modern biology laboratories. How accurate are the yours? How precise are yours? What is the difference between accuracy and precision? Find out.

Accuracy refers to performance with respect to a standard value. Precision refers to the reliability or repeatability of performance and doesn't necessarily depend on a standard at all. These two features are independent of each other. It is possible (and not so unusual) to have an instrument that is precisely inaccurate (or accurately imprecise)! Confused? Perhaps the following diagram will help.



Here, the bullseye represents the "standard" against which accuracy is judged.

To simultaneously evaluate the accuracy and precision of your pipettes you must conduct multiple measurements and compute the % Error Of The Mean and Standard Deviation of the measurements for each pipette.

A. Directions for use of Gilson Pipettes

Gilson micropipettes are virtually ubiquitous in our research labs. Often you may hear them simply referred to as "Pipettepeople".

You have 3 sizes of micropipette. The P-1000, the P-200 and the P-20. The model numbers refer to the maximum volume in microliters. Initially you will need to puzzle out which pipette to use in a given circumstance. For example, if I tell you to transfer 0.18 ml you would probably need to use the P-200 because $0.18 \text{ ml} = 180 \text{ ul}$ and because the P-200 will measure this volume with greater accuracy and precision than the P-1000 will. Generally, you chose the smallest volume pipette that will transfer the entire volume.

What would you do if I told you to measure 1.20 ml?

These pipettes are not toys. They are top of the line, expensive equipment and we demand that you use them responsibly. Here are a few simple rules.

Rules:

- Always use a disposable tip. The P-20 and P-200 use the yellow tips and the P-1000 uses blue ones. Never draw any fluid into the white barrel of the pipette itself.
 - Never lay a pipette down while there is fluid in the tip. The fluid may accidentally find its way into the barrel.
 - Never turn the adjustment scale below or above the full range settings.
 - To maximize precision, Always use the smallest volume pipette and tube for a given total volume.
-

1. Set the desired volume.

Turn the volume up just a tiny bit past the desired setting, then down.

2. Attach a tip.

Press it on firmly, with a slight twisting motion. The tip must make an air-tight seal with the pipette barrel.

3. Depress the plunger to first stop.

Most people operate the plunger with their thumb.

4. Insert the tip in the liquid you want to transfer.

Not far, just a bit below the surface.

5. Release plunger s l o w l y.

Releasing the plunger suddenly may cause the pipette to underfill.

Also, be aware that the liquid level in the tube will drop as you withdraw fluid. Make sure that it doesn't drop below the level of the tip or you will suck air.

6. As you withdraw the tip, touch it to the side wall of the tube to remove excess fluid from the exterior.


7. To dispense, depress plunger slowly to the first stop; hesitate; then depress all the way.

When pipetting viscous solutions (such as restriction enzyme stocks) it is advisable to rinse the tip by pipetting the solution up and down once or twice.

One of the most common problems for the beginner is the tendency of small liquid volumes to adhere to the pipette tip rather than dispensing into the tube. Never dispense a small volume into thin air. Always dispense into a liquid or onto the wall of a tube so that adhesion will draw the expelled liquid off the tip. Unfortunately, this may mean that the tip is contaminated and cannot be reused.

8. With the plunger still fully depressed, remove tip from the liquid.

9. Think about what to do with the pipette tip. Should you discard it or use it again? "If in doubt, throw it out."



B. Directions for checking pipette calibration.

Each member of a pair should complete the pipette calibration independently using the same pipettors.

Remember that 1.0 milliliter of distilled water weighs 1.0 gram at 25°C.

Ideally the calibration should be checked at several at points throughout its range. To save time we will check the calibration at the maximum setting only.

It is ideal to use an electronic analytical balance that is accurate to the nearest 0.1 mg. Accuracy to the nearest 1.0 mg is sufficient for checking the P-1000 and P-200, but not the P-20. (For an analysis of the way balance precision affects your calibration see the web site.)

The data and results of this calibration exercise are best recorded and presented in a TABLE. We provide these tables as an example. In subsequent exercises you will create your own tables. This is something you should do BEFORE coming to lab whenever possible.

For your P200 micropipette:

1. Use a strip of Parafilm or a plastic weighing boat on the balance pan.
2. Set the balance to zero. (The zero button may be labeled “tare”).
3. Set the pipette to its maximum volume and attach a tip.
4. Dispense the water onto the balance and record the balance reading.
5. Repeat step 4. until you have a total of five readings.
6. Compute the MEAN, the % ERROR of the MEAN and the STANDARD DEVIATION of the 5 readings. Assume 1.0 ml of water weighs 1.0 g.

The standard deviation is:

$$\text{Standard Deviation} = \sqrt{\frac{\sum [x - \bar{x}]^2}{n - 1}}$$

where “X” is each individual value, “X-bar (\bar{X})” is the average of all values, and “n” is the number of values.

The % ERROR is = $\left[\frac{(\text{MEAN VOL.} - \text{SET VOL.})}{\text{SET VOL.}} \right] \times 100$.

C. Data Sheet

P-200

Sample #	Weight Delivered (g)	Vol. Delivered (μ l)
1		
2		
3		
4		
5		

Mean = _____ μ l

% Error = _____ %

Std. Dev. = _____ μ l

D. Questions for Discussion

How much inaccuracy and imprecision are acceptable, do you think? Compare to the manufacturer's specifications given below.

The accuracy and precision of the measurements are affected by the micropipette AND by your technique. Can you devise a way to test the relative contributions of pipette error and operator error to total error?

Does the % Error Of The Mean tell you about **accuracy**, or **precision**?

Does the Standard Deviation tell you about **accuracy**, or **precision**?

Factory Performance specifications for the Pipettman-brand Autopipettes

Model	% Error of Mean	Std. Dev.
P-20	<1% @ 10-20 μ l	<0.06 μ l @ 20 μ l
P-200	<0.8% @ 60-200 μ l	<0.3 μ l @ 200 μ l
P-1000	<0.8% @ 375-1000 μ l	<1.3 μ l @ 1000 μ l

*Your pipettes may not meet these standards.

See the web site for further technical information on micropipette calibration.

E. A Test

Here are two quick tests of your ability to operate the pipettors.

- Use a P1000 to add 0.5 ml of H₂O to a microfuge tube.
Set the P200 to 100 μ l and withdraw the entire 0.5 ml.
- Use a P200 to add 0.05 ml of di H₂O to a microfuge tube.
Set the P20 to 10 μ l and withdraw the entire 0.05 ml.

III. SPECTROPHOTOMETRIC ANALYSIS OF GREEN FOOD COLOR

III.A INTRODUCTION

Each pair of students will attempt to determine the composition of green food color from a commercial set of 4 food colors (RED, BLUE GREEN and YELLOW) sold in supermarkets. "Composition" in this context means that we need to estimate the amount of each component present.

The food colors are made from pure dyes listed on the package, which are:

DYES
Tartrazine (FDC Yellow#5)
Erythrosine BS (FDCRed #3)
Allura Red AC (FDC Red #40)
Brilliant Blue FCF (FDC Blue #1)

Note that there is no GREEN DYE listed.

Your task is to determine which dye/s is/are present in green color and to estimate the molar concentration of each dye. The strategy is to determine the **absorbance spectrum** of the green food color and compare it to the known absorbance spectra of the pure dyes listed as ingredients on the package. An absorbance spectrum is a graph of the relative absorbance value (A) of the solution at different wavelengths within the visible range.

The absorbance spectra for the 4 pure dyes may be found posted on the course web site. Also see the course web site for chemical structures of the dyes and for interesting information about FDA regulation of food colors.

Please do not confuse the food *COLORS* with the pure *DYE* solutions.

III.B Procedure

- See Appendix I at the back of the Manual for a review of basic concepts in spectrophotometry, and for instructions on the use of the spectrophotometers.
- Determine Absorbance (A) for the green food color at intervals of 20 nm throughout the visible range (400 nm - 710 nm). When you find the absorbance maximum or maxima, you should take readings at 10 nm intervals on either side of the peak to define λ_{max} more precisely.

- The undiluted food color has an absorbance value > 2.0 (off scale) throughout much of the visible range. Therefore, you will need to do some systematic trial and error to find appropriate dilutions for different portions of the absorbance spectrum. The dilution appropriate for one segment of the spectrum may be different than for a different range of wavelengths. Any dilution that brings your readings below 1.5, and which is still detectable, is acceptable.
- Use the pH buffer provided to dilute the food color. The absorbance spectrum of many dyes will change with pH, so this must be controlled. Remember to mix the dilutions thoroughly before reading them. This may be conveniently done by placing a piece of "parafilm" over the mouth of the tube and then inverting several times.
- A little thought should allow you to guess the approximate (λ_{\max}) of the food color.
- Record your data in a table with four columns:

wavelength (λ)

dilution factor (DF) (for the undiluted solution DF= 1.0.)

Absorbance (A_{raw})

Corrected Absorbance ($A_{\text{corr}} = A_{\text{raw}} \times 1/\text{DF}$) [If the sample is undiluted, DF=1.0.]

- The diluted food coloring can be discarded in the drains.

III.C Results and Analysis

- Create an absorbance spectrum graph of the food color. Plot λ (nm) on the horizontal axis and Absorbance (unitless) on the vertical axis. Adjust the overall dimensions of the plot area and the scales to clearly represent the data points.

This can be done with Microsoft Excel if you choose. I can show you how.

- Determine the location of the absorbance maximum (λ_{\max}), or of the absorbance maxima if there is more than one major peak.
- Now consult the absorbance spectrum graphs of the pure dyes that we have posted on the web site. Find the dye with an absorbance maximum corresponding to each λ_{\max} value in your graph. This identifies which dye or dyes is present in your food color.

- Finally, determine the molar concentration of each dye in your food color. This requires consulting again the absorbance spectrum graphs of the pure dyes that we have posted on the web site. In these graphs, the vertical axis is the Absorbance of a 1 mM solution of the pure dye.

Knowing the absorbance of a 1 mM solution of the pure dye, the molar concentration of the dye in the food color is easily calculated on the assumption that concentration is proportional to absorbance.

Example: A 1 mM solution of Brilliant Blue FCF has an absorbance of 29 (at its λ_{max}). The Brilliant Blue FCF peak in the spectrum of a food color has an absorbance of 58 at the same wavelength. Therefore the concentration of Brilliant Blue FCF in the food color would be $58/29 = 2$ mM.

(This calculation assumes that Brilliant Blue FCF is the only component of the food color which absorbs light at this wavelength.)